ACUVUE® OASYS 1-Day

FOR ASTIGMATISM

Udobnost inspirisana suzama*[1] za **oštriji, svetliji** i **jasniji vid**[2]



Product image for illustrative purposes only.

- INKORPORIRANA VLAŽNOST[1] ZA HIDRATACIJU
- NEVIDLJIVA IVICA[1] ZA UDOBNOST
- UV ZAŠTITA[†1] ZA ZDRAVLJE OKA
- ANTIOKSIDANTI[1] PODRŽAVAJU UDOBNOST
- SOČIVO STABILIZOVANO KAPCIMA[1] ZA ASTIGMATIZAM

Šta vaši pacijenti treba da znaju o **ACUVUE® OASYS 1-Day for ASTIGMATISM**


- Zahvaljujući tome što je sočivo stabilizovano kapcima, vid je jasan i stabilan tokom celog dana, čak i pri obimnim pokretima glave i očiju.[4]
- Izuzetna udobnost[2]
- Klasa 1 UV zaštite*‡[5-7]

EYE-INSPIRED™ INOVACIJA[3]

Dizajn stabilizovan kapcima (ESD) tehnologija

ACUVUE® Contact Lenses for Astgmatism su jedina torična kontaktna sočiva sa **ESD** (dizajn stabilizovan kapcima) **tehnologijom**.[4]

- **Četiri zone stabilizacije** za prirodno poravnanje svakim treptajem za konzistentan, jasan i stabilan vid tokom dana, čak i pri obimnim pokretima glave i očiju.
- **Optička zona bez prizme**§ izbegnuta je vertikalna heteroforija koja je udružena sa simptomima kao što su vrtoglavica i mučnina.[§8,9]
- **Minimalna interakcija sa donjim kapkom** sprečava neželjenu rotaciju‖
- **Horizontalna i vertikalna simetrija** - sočiva su pod minimalnim uticajem gravitacije



For illustrative purposes only.

* Designed to support the functions of the tear film. † All ACUVUE® contact lenses have Class 1 or Class 2 UV-blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. UV transmission measured with -1.00 lens. ‡ Helps protect against transmission of harmful UV radiation to the cornea & into the eye. § Vertical heterophoria possibly caused by prism dissociation due to the presence of induced optical prism is a relevant factor for practitioners to consider when fitting toric contact lenses for monocular astigmats or those requiring a mix of toric soft contact lens designs e.g., mixed prescriptions.[8,9] Clinical studies have not been done to fully characterise the clinical effects of differences in base down prism among different contact lenses. ‖ Traditional prism and peri-ballast lens designs have a thicker area below lower eyelid than ACUVUE® Brand Contact Lenses for ASTIGMATISM that can cause interaction and may cause unwanted rotation due to lateral movement of the lower lid during blinks.

Koje su prednosti **ACUVUE® OASYS 1-Day for ASTIGMATISM**?

87% astigmatičnih očiju može se zbrinuti rasponom dioptrija brenda ACUVUE® OASYS 1-Day f or ASTIGMATISM sočiva[10].

HydraLuxe® Technology

Dizajn prožet suzama **sadrži agens za vlaženje** (polivinilpirolidon)u vidu mreže molekula sličnih suzama u visokopropusnom silikonhidrogelu koji omogućava integrisanost sočiva sa suznim filmom tokom celog dana[11,12]

Visok nivo UV zaštite

Klasa1 UV zaštite 99.9% UVB i 96% UVA[†‡5-7]

Specifikacije proizvoda

Materijal	Senofilcon A	Raspon dioptrija	**Plano to -6.00D** (0.25D steps): -0.75DC, -1.25DC, -1.75DC Axes 10-180 (10° steps) -2.25DC Axes 10-20, 70-110 and 160-180 (10° steps) **-6.50D to -9.00D** (0.50D steps): -0.75DC, -1.25DC, -1.75DC Axes 10-20, 70-110 and 160-180 (10° steps) **+0.25D to +4.00D** (0.25D steps): -0.75DC, -1.25DC, -1.75DC Axes 10-20, 70-110 and 160-180 (10° steps)
Tehnologija	HydraLuxe® Technologija Sočivo stabilizovano kapcima		
Sadržaj vode	38%	Propustljivost za kiseonik (Dk/t) [‖¶]	129x10^{-9} (-3.00D lens)
Bazna krivina	8.5mm	Markeri za orijentaciju	Vertikalne kratke linije na 6 i 12 na satu
Dijametar	14.3mm	Klasa 1 UV zaštite[†‡]	99.9% UVB, 96% UVA
Centralna debljina	0.080mm (-3.00D lens)	Karakteristike vidljivosti	Vidljiva nijansa
		Preporučeni raspored zamene sočiva	Jednokratna upotreba: zamena svakog dana

† All ACUVUE® contact lenses have Class 1 or Class 2 UV-blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. UV transmission measured with -1.00 lens. ‡ Helps protect against transmission of harmful UV radiation to the cornea & into the eye. ‖ Oxygen transmissibility at centre of a -3.00D lens using boundary-corrected, edge-corrected Dk values. Units: (cm/sec) (ml O_2/ml x mm Hg) at 35°C. Dk determined via polarographic method. ¶ Oxygen is only one of several factors related to maintaining overall eye health when wearing contact lenses.

1. JJV Data on File 2021. HydraLuxe® Technology Definition. **2.** JJV Data on File 2017. Visual Performance of ACUVUE® OASYS Brand 1-Day Contact Lenses with HydraLuxe® Technology - Consumer Use Claims. **3.** JJV Data on File 2024: ACUVUE® Brand - EYE-INSPIRED™ INNOVATIONS. **4.** Data on file 2024, ACUVUE® Brand Contact Lenses for ASTIGMATISM overall fitting success, orientation position, rotational stability and vision performance. **5.** JJV Data on File 2022. Material Properties: 1-DAY ACUVUE® MOIST, 1-DAY ACUVUE® TruEye®, ACUVUE® OASYS 1-Day with HydraLuxe® Technology and ACUVUE® OASYS MAX 1-Day with TearStable™ Technology Brand Contact Lenses. **6.** Sterner O, Aeschlimann R, Zurcher S, Osborn Lorenz K, Kakkassery J, Spencer ND, Tosatti S. Friction Measurements on Contact Lenses in Their Operating Environment: Effect of Testing Conditions on Friction. IOVS 2016 Oct;57(13):5383-5392.1. **7.** Sulley A, Mundorf M, & Arnold K. Daily Disposable Contact Lens Moduli–Are They Really Different? OVS 2016;93:E-abstract 16065 AND/OR Sulley A, Arnold K, Mundorf M. Elastic modulus of soft contact lens. Journal Japanese Contact Lens Society 2018. [Japanese]. **8.** Jackson D et al. Vertical heterophoria and susceptibility to visually induced motion sickness. Strabismus. 2012;20(1):17-23. **9.** Du Toit R et al. Tolerance to prism induced by readymade spectacles: setting and using a standard. OVS 2007;84(11):1053-9. **10.** JJV Data on file October 2023. Proportion astigmats accommodated with ACUVUE® Brand CLs for ASTIGMATISM. **11.** JJV Data on File 2018. Similarities between Mucin and Poly(N-Vinyl Pyrrolidone) (PVP). **12.** JJV Data on File 2018. HydraLuxe® Technology Definition.
2026pp15145

Sva ACUVUE® kontaktna sočiva imaju UV zaštitu kako bi pružila zaštitu od prenosa štetnog UV zračenja na rožnjaču i u oko. Kontaktna sočiva sa UV zaštitom NISU zamena za zaštitne UV apsorbujuće naočare ili sunčane naočare, jer ne štite površinu celog oka i područje oko oka. Trebali biste nastaviti s nošenjem UV apsorbujućih naočara prema preporuci zdravstvenog radnika. Prenos UV zračenja je meren na sočivu dioptrije -1.00. Za više informacija o pravilnom nošenju, nezi i sigurnosti, obratite se svom zdravstvenom radniku.

Važne sigurnosne informacije

Kontaktna sočiva ACUVUE® indikovana su za korekciju. Kao i kod nošenja bilo kog kontaktnog sočiva, mogu se razviti problemi sa očima, uključujući keratitis. Neki korisnici mogu osetiti blagu nadraženost, svrab i nelagodu. Kontaktna sočiva ne smeju se nositi u slučaju infekcije oka ili bilo kojeg drugog stanja oka, ni u slučaju sistemske bolesti koja može uticati na oko. Dalje informacije uključujući kontraindikacije, mere opreza i štetne događaje potražite u uputstvu za upotrebu.

ACUVUE®, ACUVUE® OASYS 1-Day, ACUVUE® OASYS 1 – Day for Asigmatism, HydraLuxe® su zaštićene robne marke kompanije Johnson & Johnson © Johnson & Johnson and its affiliate 2026. 2026 PP15145

Pre upotrebe pročitati uputstvo! O nameni i neželjenim reakcijama na medicinsko sredstvo posavetujte se sa lekarom ili farmacutom.

Ovlašćeni predstavnik proizvođača: Oktal Pharma doo Beograd, Vilinev ode BB, Beograd
Broj rešenja medicinskog sredstva: 001532600 2024 od 22.08.2024