ACUVUE®

Portfolio kontaktnih sočiva

EYE-INSPIRED™ INNOVATIONS[1]

Svako ACUVUE® kontaktno sočivo je inspirisano prirodnom funkcijom oka[1], na osnovu pažljivog uvida u potrebe pacijenata.

INKORPORIRANA VLAŽNOST[1] ZA HIDRATACIJU
NEVIDLJIVA IVICA[1] ZA UDOBNOST
UV ZAŠTITA[1] ZA ZDRAVLJE OKA
ANTIOKSIDANTI[1] PODRŽAVAJU UDOBNOST

In ACUVUE® Contact Lenseś

SOČIVO STABILIZOVANO KAPCIMA ZA ASTIGMATIZAM

In ACUVUE® contact lenses for Astigmatism

JEDNODNEVNA KONTAKTNA SOČIVA

NAPREDNA TEHNOLOGIJA

	SFERE	TORICI
		
	ACUVUE® OASYS 1-DAY with HydraLuxe® Technology	**ACUVUE® OASYS 1-DAY with HydraLuxe® Technology for ASTIGMATISM**
Materijal	senofilcon A	senofilcon A
Bazna krivina	8.5mm, 9.0mm	8.5mm
Dijametar	14.3mm	14.3mm
Raspon dioptrija	–0.50D to –6.00D (0.25D steps) –6.50D to –12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	Plano to -6.00D (0.25D steps): -0.75DC, -1.25DC, -1.75DC Axes 10-180 (10° steps) -2.25DC **Axes 10-20, 70-110 and 160-180 (10° steps)** -6.50D to -9.00D (0.50D steps): -0.75DC, -1.25DC, -1.75DC **Axes 10-20, 70-110 and 160-180 (10° steps)** +0.25D to +4.00D (0.25D steps): -0.75DC, -1.25DC, -1.75DC **Axes 10-20, 70-110 and 160-180 (10° steps)**
Tehnologija	HydraLuxe® Technology	HydraLuxe® Technology SOČIVO STABILIZOVANO KAPCIMA
Dk/t vrednost	121 x 10^{-9} (-3.00D lens)	129 × 10^{-9} (–3.00D lens)
UV zaštita*	Class 1, >99.9% UVB, 96% UVA[2]	Class 1, >99.9% UVB, 96% UVA[2]
Dodatne osobine	Vidljiva nijansa, „123„ oznaka pravilne strane	Vidljiva nijansa, vertikalne kratke linije na 6 i 12 na satu

VIŠEKRATNA KONTAKTNA SOČIVA

NAPREDNA TEHNOLOGIJA

	SPHERE	TORIC
		
	ACUVUE® OASYS with HYDRACLEAR® PLUS Technology	**ACUVUE® OASYS for ASTIGMATISM**
Materijal	senofilcon A	senofilcon A
Bazna krivina	8.4mm, 8.8mm	8.6mm
Dijametar	14.0mm	14.5mm
Raspon dioptrija	–0.50D to –6.00D (0.25D steps) –6.50D to –12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	Plano to -6.00D (in 0.25D steps): -0.75DC, -1.25DC, -1.75DC, -2.25DC, -2.75DC Axes 10-180 (10° steps) -6.50D to -9.00D (in 0.50D steps): -0.75DC, -1.25DC, -1.75DC, -2.25DC, -2.75DC Axes 10-180 (10° steps) +0.25D to +4.00D (in 0.25D steps): -0.75DC, -1.25DC, -1.75DC, -2.25DC, -2.75DC Axes 10-180 (10° steps) +4.25D to +6.00D (in 0.25D steps): -0.75DC, -1.25DC, -1.75DC, -2.25DC, -2.75DC Axes 10-180 (10° steps)
Tehnologija	HYDRACLEAR® PLUS Technology	HYDRACLEAR® PLUS, SOČIVO STABILIZOVANO KAPCIMA
Dk/t vrednost	147 x 10^{-9} (-3.00D lens)	129 x 10^{-9} (for -3.00D lens)
UV zaštita*	Class 1, >99.9% UVB, 96.1% UVA[4]	Class 1, >99.9% UVB, 96.1% UVA[4]
Vreme nošenja	Dnevno nošenje: zamena na 2 nedelje, Produženo nošenje: zamena na 7 dana	Dnevno nošenje: zamena na 2 nedelje, Produženo nošenje: zamena na 7 dana
Dodatne osobine	Vidljiva nijansa, „123„ oznaka pravilne strane	Vidljiva nijansa, markeri, vertikalne kratke linije na 6 i 12 na satu

All ACUVUE® contact lenses have Class 1 or Class 2 UV-blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. UV transmission measured with -1.00 lens. **
1. JJV data on file 2024: ACUVUE® Brand - EYE-INSPIRED™ INNOVATIONS. **2.** JJV Data on File 2016. ACUVUE® OASYS Brand Products Comparison. **3.** JJV Data on File 2019. Material Properties for Reusable Toric Contact Lenses, Including ACUVUE® OASYS Brand Contact Lenses for ASTIGMATISM and ACUVUE® VITA Brand Contact Lenses for ASTIGMATISM.
4. JJV Data on File 2018. Material Properties: ACUVUE® OASYS Brand Contact Lenses with HYDRACLEAR PLUS, ACUVUE® OASYS with Transitions Light Intelligent Technology, and other reusable contact lenses.
Sva ACUVUE® kontaktna sočiva imaju UV zaštitu kako bi pružila zaštitu od prenosa štetnog UV zračenja na rožnjaču i u oko. Kontaktna sočiva sa UV zaštitom NISU zamena za zaštitne UV apsorbujuće naočare ili sunčane naočare, jer ne štite površinu celog oka i područje oko oka. Trebali biste nastaviti s nošenjem UV apsorbujućih naočara prema preporuci zdravstvenog radnika. Prenos UV zračenja je meren na sočivu dioptrije -1.00. Za više informacija o pravilnom nošenju, nezi i sigurnosti, obratite se svom zdravstvenom radniku.

Važne sigurnosne informacije

Kontaktna sočiva ACUVUE® indikovana su za korekciju. Kao i kod nošenja bilo kog kontaktnog sočiva, mogu se razviti problemi sa očima, uključujući keratitis. Neki korisnici mogu osetiti blagu nadraženost, svrab i nelagodu. Kontaktna sočiva ne smeju se nositi u slučaju infekcije oka ili bilo kojeg drugog stanja oka, ni u slučaju sistemske bolesti koja može uticati na oko. Dalje informacije uključujući kontraindikacije, mere opreza i štetne događaje potražite u uputstvu za upotrebu.
ACUVUE®, ACUVUE ® OASYS 1-Day, ACUVUE® OASYS 1-Day for Astigmatism, ACUVUE® OASYS for Astigmatism, HydraLuxe®, Hydraclear®PLUS su zaštićene robne marke kompanije Johnson & Johnson © Johnson & Johnson and its affiliate 2026. 2026 PP14883

Pre upotrebe pročitati uputstvo! O nameni I neželjenim reakcijama na medicinsko sredstvo posavetujte se sa lekarom ili farmaceutom

Ovlašćeni predstavnik proizvođača: Oktal Pharma doo Beograd, Vilinev ode BB, Beograd
Broj rešenja medicinskog sredstva: 001532600 2024 od 22.08.2024